



Refrigerated Centrifuge

INV- RC15



Description & Features

- As one of the pursuit of ESG values, it supports energy saving functions.
- Standby Off function, electric lid-locking system
- Eco Shut-off setting for refrigeration
- Unnecessary energy use can be avoided.
- Solitary port for water collection in chamber
- Through efficient refrigeration control technology and internal design that minimizes the frictional heat of the rotor, it guarantees safe testing of temperature -sensitive samples.
- User-friendly 4.3" full color touch LCD controller helps to set the parameter at modern way. All set points are visible at a glance and setting is easy. - Touch pad lock doesn't allow misuse during operation avoided.
- Emergency lid lock release in case of power failure to take samples out.
- Unauthorized access is prohibited through password
- All the rotor are metallic and Autoclavable

INNOVA BIO LIFE TECHNOLOGY
1st Floor, H. No. 13/10, Ram Nagar
Karnal, Haryana-132001





Specifications

Max RPM	15,000RPM(Angle Rotor)	5,000RPM(Swing Rotor)
Max RCF	23,847 xg	4,947 xg
Max capacity	85ml x 6	750ml x 4
Accel / Decel phase	0 ~ 9 (10set)	
Set Temp.	-10°C ~ +40°C (1°C increments)	
Ref. System	Non-CFC	
Noise Level	≤60dB	
Set Time	99 hour 59 min 59sec or continuous	
Time Count Program	From start / From set-up RPM	
Program	100 memory	
Rotor ID	Automatic Recognition	
Display	4.3" LCD Touch Screen	
Display message on screen	RPM, RCF, Time, Door open/close, Temp ACC, DEC, Rotor ID, Program, Comp. On/Off, key lock	
RPM Range	400 to 15,000 rpm	
Power supply	AC 230V ± 10%, 50-60Hz, 1Phase	
Power	2,5KW	
Dimension	765 x 687 x 375 mm (H:915mm with open lid)	
Weight	95kg	

Available Rotors

Rotor no.	Capacity	RPM
A480.2T	Fixed Angle Rotor48 x 0.2ml	12,000
A481.5T	Fixed Angle Rotor 48 x 1.5ml	13,000
A301.5T	Fixed Angle Rotor30 x 1.5ml	15,000
AS301.5TX	Fixed Angle Rotor30 x 1.5ml (aerosol-tight)	15,000
A1215T	Fixed Angle Rotor 12 x 15ml	12,000
A1215TC	Fixed Angle Rotor12 x15ml(conical)	13,000
A650T	Fixed Angle Rotor 6 x 50ml	12,000
A650TC	Fixed Angle Rotor 6 x 50ml(conical)	13,000
A850T	Fixed Angle Rotor 8x 50ml	12,000
A685T	Fixed Angle Rotor 6 x 85ml	12,000
S4250TE	Swing Out3mlx12/48	4000
S4250T	Swing Out5mlx8/32	4,500
S4750T	Swing Out 15ml Conicalx14/56	4000
	Swing out 50mlx5/20	4000
S4100T	Swing out1.5mlx6/24	5000
S47501	Drum rotor60x1.5/2ml or 120x0.4ml	14000

INNOVA BIO LIFE TECHNOLOGY
1st Floor, H. No. 13/10, Ram Nagar
Karnal, Haryana-132001

